



TECHNICAL SHEET



BSU154 OB2012 FA WHITE MEDICAL SUCTION UNIT WITH WALLBRACKET

DESCRIZIONE

- Medical Suction Unit OB 2012 powerful, easy and compact.
- Silent and powerful double piston motor of more than 30 lt./min.
- Adjustable linear suction, with control-manometer from 0 to > 800 mbar.
- Autonomy of about 45 minutes.
- New box made in ABS/UL.
- Autoclavable 1000 ml-jar for the secretion collection with overflow valve and filter directly integrated in the cover.
- Patient tube with Jankauer.
- Cable for the connection to the 12 Volts C.C. of the vehicle.
- Dimensions: about 35x12x24 cm.
- Weight: about 4,5 Kg.
- Manufactured in compliance with ISO 10079-1:1999, EN 1789:1999 norms and classified in class IIb.
- CE certified under the MDD93/42/EEC, "e1" 23555
- homologated according to 95/54/EC directive.
- Wall bracket anchorage which integrates the recharging system, tested at 10G in compliance with EN 1789:1999 regulation.



CLASSIFICATION: The device is all in class IIb accordingly to MDD 93/42/CEE

QUALITY CONTROLL: 100% of the entire supplying

PACKAGING: Singularly packaging

DEVICE CONFORMITY / REFERENCE NORMS

MDD47/2007/EC	Directive medical devices
ISO14971	Application for risk managing on medical devices
EN980	Graphyc symbols used for the labelling of medical devices
ISO10079-1	Vacuum cleaners for medical use. electronically powered suction. Safety requirement
EN1789	Vehicles and their medical equipment. Ambulances

DOCUMENT REFERENCE		
ADITED	APPROVED	DATE
DI/RAQ PERINI LAURA	DI/AD BRAZZO DANIELE	07/05/2010



Via E. Ferrari, 29 – 39100 Bolzano Italy
Tel. +390471932893 – Fax +390471932331
www.boscarol.it – Email info@boscarol.it