



5 OB SSPLINT SET

- Splint system set with aluminium core for easy limbs immobilization in case of fracture or sprain
- Made in neoprene material
- Adaptable to fit each body frame
- It does not cause tissue compressions
- It maintains body temperature
- It reduces congestion around injured areas
- Velcro modelling closing
- Lightweight and easily washable even in washing-machine
- Least encumbrance
- Complete with bag
- Dimensions with bag: about 64x15x29(h) cm
- Weight: about 1,8 kg
- Not suitable for femur fractures
- Ref: IMM124030

SSPLINT SYSTEM SET WITH ALUMINIUM CORE
FOR ESAY LIMBS IMMOBILIZATION



LIGHTWEIGHT AND EASILY WASHABLE
EVEN IN WASHING-MACHINE

