



TRANSFER SHEET

- New transfer sheet, even more resistant
- Made of red tearproof impermeable PVC material
- 8 side handles
- Easily washable and disinfectable
- Dimensions: about 187x70 cm
- Weight: about 1200 g
- Loading: about 150 kg
- In accordance to EN1865 norm.
- Perimeter reinforced with a 4cm, double sewn band on both sides
- In accordance to EN597-1 on devices flammability
- Ref: IMM120321



PERIMETER REINFORCED
WITH A 4CM,
DOUBLE SEWN BAND
ON BOTH SIDES



EVEN MORE RESISTANT

