



Learn more about our new OB 3000

MODEL	OB 3000	OB2012
Flow rate	33 LPM	30 LPM
Max vacuum	800 mbar (80 kPa)	800 mbar (80 kPa)
Rated voltage	11÷25 Vdc	12 Vdc
Autonomy	> 70 minutes	> 60 minutes
Weight	2,6 Kg	4,2 Kg
Dimensions	360 x 244 (h) x 195 mm	370 x 290(h) x 130 mm
IP rating	IP44	IP34d
Battery type	removable LIPO 5 Ah 11,1 V	12 V 4 Ah Pb/Ca (lead-acid battery)
Charging time	80% charge: 6 hours 100% charge: 10÷15 hours	100% charge: 15 hours
Warranty	5 years	3 years
Periodic maintenance	2 years (automatic reminder)	2 years
IEC60601-1-12	Approved	Approved
RTCA DO160G	YES – AVIO version	NO



Visit our website to learn more about special launch prices and promotions.

Please consult the user manuals of suction devices in order to have complete technical data.



Learn more about our new OB MINI

MODEL	OB MINI	OB1000
Flow rate	24 LPM	24 LPM
Max vacuum	850 mbar (85 kPa)	850 mbar (85 kPa)
Rated voltage	11÷25 Vdc	12 Vdc
Autonomy	> 100 minutes	> 45 minutes
Weight	1,8 Kg	3,3 Kg
Dimensions	290x 240 (h) x 95 mm	335 x 290(h) x 105 mm
IP rating	IP43	IP34d
Battery type	LIPO 5 Ah 11,1 V	12 V 3.4 Ah Pb/Ca (lead-acid battery)
Charging time	80% charge: 6 hours 100% charge: 10÷15 hours	100% charge: 15 hours
Warranty	3 years	2 years
Periodic maintenance	2 years (automatic reminder)	2 years
IEC60601-1-12	Approved	Approved
RTCA DO160G	YES – AVIO version	YES – AVIO version



Visit our website to learn more about special launch prices and promotions.

Please consult the user manuals of suction devices in order to have complete technical data.

