

# OB 3000



## NEW Medical Suction Unit OB 3000

Max. suction rate 800 mbar (80 kPa)

Max. suction flow 33 liter/minute

Weight < 2,8 kg

1 lt LINER version with disposable bag

and FA version with autoclavable jar

Lithium polymer battery (solid-non organic)

Battery status indicator with no need to switch on the device

Battery to be replaced by user

Works with 12 and 24 Vcc



Inner microprocessor for all data record and functioning tracing  
Safety inspection reminder

Analogic fluorescent vacuum meter, with mmHg scale upon request  
OB 30WB wall bracket (EN1789 conform) and universal main supply  
(from 90 to 240 Vac)

Compliant to all current norms, 93/42/CEE Directive and main reference  
norm ISO10079-1



# SPECIAL LAUNCH PRICE

valid until May, 31st 2016

~~979,00 €~~

**490,00 €**

All prices are net, taxes excluded and without further discount.

# OB 3000

## New Medical Suction Unit OB 3000

Ref. and technical data



Features	OB3000
Aspiration flow rate***	33 LPM FAF
Max vacuum*	800 mbar ±10%
Power source	11÷25 Vdc
Internal battery	YES
Battery type**	LI-PO
Extractable battery	YES
Operating Autonomy	> 60 minutes
Vacuum meter****	<ul style="list-style-type: none"> <li>• Analogic</li> <li>• Double scale</li> <li>• Fluorescent</li> </ul>
Recharging bracket	OB20WB
Vacuum regulation	MECHANICAL
Autonomy indicator	<ul style="list-style-type: none"> <li>• 4 leds with microprocessor</li> <li>• Autonomy test</li> </ul>
Secretion jar	FA (1000 cc) or LINER (1000 cc)
Internal maintenance and data managing	YES Automatic visual alarm (LED)
Transport bag	NO
Shoulder strap	YES (integrated in the device)
Dimensions	340x110x220 mm
Weight	<2,8 kg
Warranty	5 years



\* The atmospheric pressure affects maximum vacuum rate  
 \*\* PB/CA lead/calcium battery – LI-PO lithium polymer battery (certification UN38.3 – IEC62133), suitable for air/sea transport  
 \*\*\* FAF = free air flow (reference norm ISO10079-1)  
 \*\*\*\* Available with mmHg scale upon request

<b>BSU3000</b>	OB 3000 FA MEDICAL SUCTION UNIT
<b>BSU3020</b>	OB 3000 FM MEDICAL SUCTION UNIT
<b>BSU800</b>	OB2012 WALL BRACKET (Compatible with OB 3000)
<b>BSU892</b>	POWER SUPPLY 2P - EU

