



**VAL43141**

## **BOSCAROL HIGH IMPACT CASE M**

- Boscarol High Impact small orange case with high resistant thick walls
- Automatic pressure release valve
- Easy to open double throw latches
- Load tested soft grip handle
- Holes for double padlocking
- Holes for shoulder strap
- Stackable interlocking system
- Thick body material
- Extra fixing points
- Fully length hinges
- Nylon pins, corrosion free
- Pick N Pluck foams in the bottom and convoluted foam available on demand
- Screen printing space on the lid
- Rectangular resin label space on the lid
- Dimensions internal: about 42,6x29x15,9 cm
- Dimensions external: about 46,4x36,6x17,6 cm
- Weight: 2,78 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.