



VAL43142

BOSCAROL HIGH IMPACT CASE L

- Boscarol High Impact small orange case with high resistant thick walls
- Automatic pressure release valve
- Easy to open double throw latches
- Load tested soft grip handle
- Holes for double padlocking
- Holes for shoulder strap
- Stackable interlocking system
- Thick body material
- Extra fixing points
- Fully length hinges
- Nylon pins, corrosion free
- Pick N Pluck foams in the bottom and convoluted foam available on demand
- Screen printing space on the lid
- Rectangular resin label space on the lid
- Dimensions internal: about 50x35x19,4 cm
- Dimensions external: about 55,5x42,8x21,1 cm
- Weight: 3,86 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.