



**BOR87010N**

## **EMERGENCY RUCKSACK MAXIMED 2.0**

- Maximed Rucksack 2019 Edition made of red nylon
- Tearproof and waterproof material, allowing the use of the rucksack during rainy conditions
- Features a reflective stripe on the padded shoulder straps for improved visibility, as well as on the front and side pockets
- The part in contact with the back is equipped with a shock absorbing cushions
- The rucksack is equipped with two front pockets
- 1 empty external large pocket and 1 external large pocket equipped with 4 removable coloured pouches
- The 4 coloured inner pouches have a transparent front and are removable
- Special straps on the inside for secure fixing of an oxygen cylinder
- Dimensions 48x30x52cm
- Capacity: about 54 litres
- Weight: about 4.5kg
- To be hand washed at a maximum water temperature of 30°C
- Provided empty, but with the possibility to be equipped with medical devices on request
- The Maximed Rucksack complies with the REACH regulation and it is not a medical device

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.