



**BOR81146N**

## **EMERGENCY RUCKSACK ATLS MODULE 2.0**

- ATLS modular rescue backpack
- Made of tear-proof and waterproof Nylon
- Reflective strip on padded shoulder straps for better visibility
- Transparent front window for logos
- Large front pocket allows safe storage of large medical devices
- Soft back pads to reduce back fatigue
- Chest strap for easier carrying
- Rigid shockproof polypropylene tub with 4 rigid subdivisions for storing bags
- 4 inner removable coloured pouches with transparent fronts
- Dimensions: approx. 32x25x53 cm
- Volume: approx. 35 l
- Weight: approx. 3.3 kg
- Maximum loadable weight 30 kg
- Sold empty but can be equipped according to specific customer requests

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.