



**BOR87050**

## **RESPONDER BAG RESCUE BACKPACK**

- RESPONDER BAG rescue backpack
- Tear-proof and waterproof red fabric
- Equipped with:
  - top handle for easy carrying
  - large, full-opening front pocket with 1 strip of elastic on the inside and 2 on the outside for securing materials
  - 2 small side pockets
  - front and side reflective strips
  - 6 removable coloured pouches with transparent front and Velcro-type strips
  - reinforced bottom and rubber feet
- The part in contact with the back is equipped with shock-absorbing cushions
- Supplied empty but can be fitted with medical devices at customer's request
- Can be hand-washed at a maximum temperature of 30° C
- Dimensions: approx. 42x23x46 cm
- Weight: approx. 1.5 kg
- Maximum transportable weight 20 kg

DIR/QM – Rev. 01\_26/07/2023

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.