



**BOR87230**

## **RUCKSACK MAXI TROLLEY RED EMPTY**

- Ideal for transporting material in hospital wards
- Tear-proof and waterproof red PVC fabric
- Equipped with:
  - 2 reflective strips
  - retractable trolley handle
  - retractable padded shoulder straps
  - 2 equipped front pockets
  - 3 side pockets
  - 8 removable coloured pouches with transparent front
- Dimensions: approx. 52x29x61 cm
- Weight: approx. 5.2 kg
- Supplied empty

*DIR/QM – Rev. 01\_16/08/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.