



**126140107191**

### **YANKAUER RIGID SUCTION PROBE COMPLETE\***

- Yankauer rigid suction probe, complete of 180 cm long connecting tube with female connection
- Patient aspiration using only one hand
- Rigid cannula instantly ready and easy to use
- Vacuum regulation by the thumb of the hand
- Suitable for every kind of medical secretion suction unit
- Singularly sterile packaged
- EC medical device

*DIR/QM – Rev. 02\_18/07/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.