



**BSU504**

## **SECRETION JAR 500 ML**

- Suction jar for medical suction units, made in autoclavable and shockproof Macrolon with pressure closing lid made in silicone-rubber
- Complete with overflow valve
- Capacity: 500 ml

*DIR/QM – Rev. 02\_12/06/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.



For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.

**OSCAR BOSCAROL S.r.l.**  
Via Enzo Ferrari, 29 - 39100 Bolzano BZ Italy  
VAT Number e Fisc. Code 01458460217  
Ph. +39 0471932893 Fax. +39 04711880333  
info@boscarol.it - www.boscarol.it

Printed: 03/04/2025 04:31:20