



4501051

## PRESSURE REDUCER P300 UNI

- Pressure reducer for medical oxygen compliant with the standard ISO 10524-1
- Made of surface treated brass
- P1 (maximum inlet pressure): 300 bar
- P2 (nominal output pressure): 4,5 bar
- Reducer capacity at 200 bar (20.000 kPa): max. 250 lpm (litres per minute)
- Oxygen cylinder connection type UNI 4406
- Output connection: 3/8" (three eighths of an inch)
- Compliance: CE marking according to MDD 93/42/CEE
- Oxygen compatible according to ISO 15001
- Operating temperature: -20÷60 °C
- Storage and conservation temperature: -20÷60 °C
- Warranty: 2 years from date of purchase (see warranty conditions)
- Periodic safety check: 5 years from the date of manufacture
- Dimensions: about 110x51x58mm
- Weight: about 580g
- Manufacturer: Hersill S.L. – Spain
- Distributed by Oscar Boscarol Srl, Via Enzo Ferrari 29 – 39100 Bolzano, Italy

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.