



CH 30 - 10 mm

185420-000300

### WENDL'S CANNULA CH 30 10MM

- Wendl nasopharyngeal cannula, made of WIRUPREN<sup>®</sup> with movable disc and enlarged tip
- For nasopharyngeal intubation in anesthesia, intensive care and emergency.
- Supplied individually in sterile packaging
- Size CH 30
- Inner diameter: 7.0 mm
- Outer diameter: 10.00 mm
- Length: 170 mm
- CND CODE: R010101
- RDM CODE: 627420

*DIR/QM - Rev. 01\_06/07/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.