



**RIA50822**

### **INTUBATION BAG OPTICAL FIBER\***

- Intubation bag complete of:
- 1 fiber optic laryngoscope with 3 blades (1-2-3)
- 1 Magill forceps
- 1 Pean forceps
- 1 scissors
- 5 endotracheal tubes
- 3 Guedel airways
- 1 syringe
- 1 mouth opener
- 1 plaster roll
- 2 Wendl airways
- 1 stylet
- Dimensions: about 32x22x7,5(h) cm

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.