



**RIA53204**

#### **LARYNGO. O.F. MC BLADES SIZE 1/2/3/4**

- Mc Intosh fiber optic laryngoscope
- Made in stainless steel
- 2,5V halogen lamp set inside the blade for an easy sterilization
- Accurately finished off
- 4 Mc Intosh F.O. blades (1-2-3-4), standard handle
- Supplied with a practical roll case
- The handle is designed to contain 2 C batteries (not included)

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.