



394600

### 3-WAY TAP

- Luer-Lock fittings, 1 male swivel outlet and 2 female inlet, hermetically sealed up to 44Psi (300kPa)
- Faucet outlets protected by Luer-Lock caps until used.
- Reduced length female Luer-Lock fittings to nullify debris during discontinuation of infusion therapy.
- Anti-turbulence Luer-Lock fittings to reduce the risk of microbubble formation inside the tap.
- Resistance to rupture within the limits of use and within existing guidelines.
- Possibility of connection to two or more infusion systems.
- TYVEK packaging.
- Sterilized by irradiation in accordance with existing regulations.
- Apyrogenic, disposable, individually packaged.
- Non-toxic materials compatible with pharmaceutical products (drugs, disinfectants, infusion solutions) with which they may come in contact during normal use practice
- Storage conditions: at room temperature, in a cool and dry place away from direct light and heat sources.
- Expiration date: 3 years from the production lot date printed on the package.
- Resistant to lipids, blood and blood products.
- Product stability over time: within the validity limits indicated on the package
- **EXPIRY 12/2024**

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.