



## ETAN30C02

### 3.0MM CUFFED ORO-NASAL TRACHEAL TUBE

- Tracheal, oro-nasal tubes, cuffed, made of polyvinyl chloride (PVC), thermosensitive, intended for the following uses:
- Positive pressure ventilation;
- Induction of anaesthesia when it is not possible to use a facial mask;
- Removal/protection of secretions from the trachea and bronchi during aspiration. Single-use. Oro-nasal;
- Radiopaque line along the entire length;
- 15M connector;
- Centimetre scale;
- Low-pressure, high-volume atraumatic cuff;
- Blue indicator balloon with Luer connector and unidirectional valve;
- Murphy eye with rounded, atraumatic edge;
- Depth markings;
- Atraumatic, pre-curved Magill-type tip in accordance with applicable ISO/EN standards.

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.