



**PF-273P-0916**

## **HUMIDIFIER FOR FLOWMETER INPUT 9/16"**

- Humidifier for oxygen therapy with 9/16" inlet
- Jar and cap autoclavable at 121° C - 2bar - max. 20 minutes
- Silicone inner tube with terminal equipped with a filter for homogeneous diffusion of the gas in the water
- Hose connector positioned on the humidifier cap for connection to patient-supply-devices (masks, nasal catheters...)
- Indication of minimum and maximum water filling limits
- Latex safe/free
- Made in Italy
- Class IIa medical device in compliance with 93/42/EEC as amended and UNI 11341 standard.

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.