



32693

## YTON ANEROID SPHYGMOMANOMETER WITH PHONENDO

- KIT YTON consisting of sphygmomanometer and stethoscope
- Classic sphygmomanometer with nylon cuff, Velcro fastener, latex pear with chrome valve and detachable manometer which can be attached to the cuff.
- Chromium-plated air valve with air release control via screw valve.
- Scale graduated every 2 mmHg, from 0 to 300 mmHg.
- Manual control of air release possible from 260 mmHg to 15 mmHg.
- Deflation time max. 10 sec.
- Floating zero (pointer does not have a stop point but oscillates freely)
- Supplied with a soft leatherette bag with zip, in a cardboard box.
- Colour: black
- Class Im medical device in compliance with Directive 93/42/EEC and subsequent amendments and supplements

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.