



**DIA100070**

## **PEDIATRIC BIG PALMAR SPHYGMOMANOMETER**

- Aneroid palmar sphygmomanometer
- Big manometer with graduated scale up to 300 mm/Hg
- Different colours manometer: red and blue; red for an immediate recognition of possible hypo- or hypertension conditions
- Complete with 2 cuffs, newborn and paediatric in fantasy fabric material with velcro closing
- One-handed operations allowed, accurate regulation valve
- Supplied with fantasy fabric case
- Manometer Ø 70 mm

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.