



## RM300

### DIGITAL THERMOMETER

- Measurement of body temperature by oral, axillary or rectal contact
- Measurement in 60 seconds
- Units of measurement in degrees Celsius (°C) and Fahrenheit (°F)
- Sensed temperature range: 32°C to 42.9°C
- LCD display with 3 displayed digits
- Water resistant tip (IP27 degree of protection)
- 1.5 V alkaline battery. Autonomy of approx. 200 hours
- Keeps track of the last measurement
- Class IIa EC medical device in compliance with MDD/93/42/EEC and and subsequent amendments and additions.

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.