



HS-30Q

PROFESSIONAL CHILDREN STETHOSCOPE

- Professional stethoscope ideal for use in the paediatric sector
- High sensitivity and excellent sound transmission
- The head of the device is made of chromed metal
- Ear tube made of durable latex-free plastic material, easily cleanable and disinfected
- CE Class I medical device according to European Regulation 2017/745 (MDR)
- Imported and distributed by Oscar Boscarol SRL

DIR/CEO – Rev. 02_04/11/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.