



HS-50A

ANEROID SPHYGMO WITH STETHOSCOPE

- Kit consisting of sphygmomanometer and stethoscope
- Classic sphygmomanometer with nylon cuff, Velcro fastener, latex pear with chrome valve and detachable manometer which can be attached to the cuff.
- Chromium-plated air valve with air release control via screw valve.
- Scale graduated every 2 mmHg, from 0 to 300 mmHg.
- Manual control of air release possible from 260 mmHg to 15 mmHg.
- Floating zero (pointer does not have a stop point but oscillates freely)
- Supplied with a soft leatherette bag with zip, in a cardboard box.
- Colour: black
- EC medical device

DIR/CEO – Rev. 02_04/11/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.