



## HS-201X

### SPHIGMOMANOMETER PALMAR WHITE

- Palm-sized sphygmomanometer made of white ABS
- Very light medical device that can be easily used by everyone (approximately 250 g)
- Large dial with graduated scale up to 300 mm/Hg
- Nylon bracelet with Velcro closure for quick application
- Measurement accuracy:  $\pm 3$  mmHg on the measurement taken
- Practical protective and transport case
- CE medical device in class Im
- Imported and distributed by Oscar Boscarol SRL

DIR/CEO – Rev. 02\_04/11/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.