



**HS-GF-301**

## **PALM SPHYGMOMANOMETER WITH 3 CUFF SET**

- Professional handy palmar sphygmomanometer
- Wide manometer with graduated scale up to 300 mm/Hg
- Complete with 3 cuffs
- Pediatric cuff range 18,4-26,7cm
- Adult cuff range 25,4-40,6cm
- Obese cuff range 33,3-51cm
- Supplied with case

*DIR/CEO – Rev. 02\_04/11/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.