



**BOR84065**

### **E-HOLSTER RED COMPLETE\***

- Complete belt E-Holster
- Made in tearproof and water-resistant red material
- 5 compartments with pockets for gloves and for Pen Light
- Belt loop of about 5,5 cm
- Dimensions: about 14x3x19(h) cm
- Complete with:
  - 1 first aid scissors
  - 1 anatomic forceps
  - Cocher forceps
  - 1 Pen Light
  - 1 Tourniquet

For further information and/or technical data related to the product, please refer to the Oscar Boscrol Company.