



**FCT2000**

## **STANDARD DISPOSABLE ICE 13.5x17 CM**

- Disposable ice with controlled temperature (0°C for about 30 minutes)
- Its immediate use during traumas avoids haematoma and swelling formation
- It stops little haemorrhages and considerably reduces the pain
- It can be kept without limits at ambient temperature
- Made in PVC material
- Dimensions cm 17x13,5
- CE certified in risk class IIa
- CND code M9001 – RdM code 163045/R

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.