



**VAL44060**

## **PERIPHERAL VEINS CANAL SET**

- Peripheral veins cannulation set in pouches with transparent frontal, containing:
- 1 haemostatic belt
- 5 canula needles (16/17/18/20/22)
- 3 butterfly needles (21/23/25)
- 1 silky plaster roll 2,5 cm
- 4 disinfectant wipes
- 2 IV administration sets
- Dimensions : about 27x12x9(h) cm

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.