



## LHIOD1005

### LH IODO 10 500ML DISINFECTANT

- Broad-spectrum disinfectant for skin disinfection and hygiene in wards at risk of infection.
- Antiseptic-antibacterial for use in the following cases:
  - Preparation of the operating field
  - Skin disinfection prior to injection therapy
  - Use as an antiseptic-microbicide in medical surgeries.
- Surgical medical device

DIR/QM – Rev. 01\_04/06/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.