



9015PED

FAST COLLAR - PEDIATRIC

- The Fast Collar is the new cervical collar adjustable to different heights.
- Allows immobilization of the cervical spine and minimizes flexion-extension.
- Ideal for those who want to have, with a single article, several possibilities of immobilization.
- Ideal for post-traumatic situations in the cervical area, where a valid support and/or protection of the skeletal system and the muscle bands is necessary.
- Central hole for the control of the trachea and the execution of a tracheotomy procedure.
- Has been developed for use during medical emergency and flat one-piece design allows efficient placement in limited spaces.
- Polyethylene structure, black LD padding, nylon buttons and Velcro strap latch.
- If used when skin lesions occurred, protect the wound with an appropriate bandage; in such situations, do not use the device in direct contact with the wound.
- Radiotransparent - Latex-Free - Unsterile - Single-use - Disposable
- The device is disposable. However, it is possible to wash in water at 30°C with neutral soap or detergent.
- RDM Code 611391/R
- CND Code Y060312

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.