



**IMM120180**

## **SET OF 2 VACUUM SPLINTS+ MINIPUMP+BAG**

- Set of 2 vacuum splints, leg and arm
- Suitable for paediatric use
- Absolute, fast and painless immobilisation of fractured limbs (even exposed ones)
- High resistant washable material
- Compact suction plastic valve with external  $\frac{1}{4}$  on/off.
- Ensuring a perfect immobilisation but do not compress, and allow physiological blood circulation
- 100% X-ray transparent (no metal parts included)
- Complete with transport bag and mini-pump
- Arm splint dimensions: about 47x75 cm
- Leg splint dimensions: about 64x92 cm
- Weight: 3,6 Kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.