



**IMM120184**

### **3 VACUUM SPLINTS SET+ PUMP+BAG**

- Set of 3 vacuum splints: hand, leg and arm
- Suitable for paediatric use
- Absolute, fast and painless immobilisation of fractured limbs (even exposed ones)
- High resistant washable material
- Compact suction plastic valve with external ¼ on/off.
- Ensuring a perfect immobilisation but do not compress, and allow physiological blood circulation
- 100% X-ray transparent (no metal parts included)
- Complete with transport bag and mini-pump
- Hand splint dimensions: about 50,5x33 cm
- Arm splint dimensions: about 47x75 cm
- Leg splint dimensions: about 64x92 cm
- Dimensioni of set: about 54x24x44(h) cm
- Weight: about 4 Kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.