



**IMM123120**

## **2 RIGID SPLINTS SET WITH BAG**

- 2 rigid splints set for upper and inferior limbs immobilization
- Made in PVC and Nylon
- Washable and X-ray transparent
- Velcro closing
- Transport bag with handle
- Dimensions inferior limb: about 74x40x3,5 cm
- Dimensioni upper limb: about 51x27,5x3,5 cm
- Dimensions in bag: about 76x19x19 cm
- Weight: about 2,1 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.