



**IMM124030**

### **SET OF 5 OB SSPLINTS WITH BAG**

- Splint system set with aluminium core for easy limbs immobilization in case of fracture or sprain
- Made in neoprene material
- Adaptable to fit each body frame
- It does not cause tissue compressions
- It maintains body temperature
- It reduces congestion around injured areas
- Velcro modelling closing
- Lightweight and easily washable even in washing-machine
- Least encumbrance
- Complete with bag
- Dimensions with bag: about 64x15x29(h) cm
- Weight: about 1,8 kg
- Not suitable for femur fractures

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.