



**IMM120214**

## **ARM VACUUM SPLINT BASIC+TRANSPORT BAG**

- Anatomic conformed vacuum splint suitable for leg and ankle
- Absolute, fast and painless immobilisation of fractured limbs (even exposed ones)
- High resistant, easily washable material
- Compact suction valve with external ¼ on/off
- Protection bags in tearproof, impermeable material, with rapid hook/unhook fixing and moulding system
- Ensuring a perfect immobilisation but do not compress, and allow physiological blood circulation
- 100% X-ray transparent (no metal parts included)
- Dimensions: about 68x38x53 cm
- Weight: about 0,8 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.