



## HS201C1

### HANDHELD SPHYGMOMANOMETER\*

- Hand-held sphygmomanometer with black cuff
- Large pressure gauge with graduated scale up to 300 mm/Hg
- Nylon cuff with Velcro closing for rapid use
- Inflation bulb with pressure regulation valve
- Measuring accuracy:  $\pm 3$  mmHg
- Supplied with transport case
- Risk class Im, according to MDD 93/42/CEE (measuring instrument)

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.